


Some of our glass laboratory electrodes have gotten an update! Use the table below to identify replacements for part numbers you recognize and check out new additional features.

Part Numbers		Features			
New Part Number	Replaces	Junctions	Cap Type	Strain Relief	BNC Boot
pH1400	SG201C	Single	Standard	No	No
pH2400	S200CD SG200C	Double	Sensorex Freedom Cap	No	Yes
pH2500	SG350CD SG350C	Double	Sensorex Freedom Cap	No	Yes
pH3400	SG1041CD	Double	Refillable Pump	Yes	Yes



= *New Features*