




PLASTIC LABORATORY pH ELECTRODES PART NUMBER EQUIVALENTS GUIDE

Some of our plastic laboratory electrodes have gotten an update! Use the table below to identify replacements for part numbers you recognize and check out new additional features.

Part Numbers		Features				
New Part Number	Replaces	Body Material	Junctions	Cap Type	Strain Relief	BNC Boot
pH1000	S150C	Polycarbonate	Single (1)	Standard	No	No
pH2000	S200CD	Ultem®	Double (2)	Sensorex Freedom Cap	No	Yes
pH2100	S350CD	Ultem®	Double (2)	Sensorex Freedom Cap	No	Yes
pH2200	S175CD	Ultem®	Double (2)	Sensorex Freedom Cap	No	Yes
pH3000	S1021CD	Ultem®	Double (2)	Refillable Pump	Yes	Yes



= *New Features*